



Stepping into Room Temperature Hydrolysis for a Comprehensive Panel

Mauricio Vargas, PhD
IACFT 2021

Outline

1. About Kura and Finden
2. What is B-One?
3. How does B-One perform under:
 - Optimum conditions?
 - Challenging conditions?
4. Conclusions

From Chilean Patagonia to the world



ISO 9001:2015
certified



Present in +15
countries



State of the art
R&D and
production lab in
Chile



Offices in the
United States



Providing the
world largest
laboratories

**+25 million
people tested
with KURA
Biotech
reagents
each year.**


KURA's Way

Inspired by nature. Moved by science.

We manufacture our enzyme in a very unique way by creating **our own genomic database from nature.**

KURA is using the most advanced tools afforded by **modern biotechnology and computer-modeling, modifying enzymes through directed evolution** to achieve more specific results in the most demanding analysis.

The perfect fusion between nature, scientific approach and expert knowledge.



KURA
biotech

Our Brands

KURA
biotech

TOXICOLOGY

 **finden**
KURA


GENETICS

a-venire
KURA

PROTEOMICS


blikka
KURA

KURA
biotech



Let's talk about Finden and β -Glucuronidase hydrolysis

Classic β -Gluc Hydrolysis




What is B-One?

B-One® High-efficiency recombinant β -Glucuronidase for RT Hydrolysis



- **>85% Recovery**
- **Codeine-6- β -D-Glucuronide (2,500 ng/mL)**
- **5 minute incubation**
- **At Room Temperature**
- **Stable to Store at Room Temperature**

- All-in-One β -Glucuronidase stabilized in Buffer Solution
- Highly purified, clean enzyme, Dilute & Shoot friendly
- No additional reagent mixing or clean up is needed
- Allows fully automated sample preparation
- Reduced labor and fewer processing errors thanks to simplified protocol




B-One® is highly stable

Stability measurements:

- Up to 3 months if left **on the bench**
- Up to 1 year if stored **refrigerated**

Eliminate storage concerns since it is stable at room temperature and refrigerated.





B-One's challenge:

performance in real world conditions

What if there are
other conjugated
drugs in the
sample?

I usually have
higher analyte
concentrations.



What if I
change my
extraction
method?

What if I take
longer than 5
minutes?

Goal

Test B-One
under real UDT
conditions



We listened,
transferred the
information,
and applied it.


New Conditions to be Tested






Results

1st A Simple Experiment to Start



Recovery under a simple protocol



EASY TO CLEAVE


Temazepam Glucuronide

Amitriptyline-N-β-D-Glucuronide

(+)-11-nor-9-Carboxy- D9-THC glucuronide

Oxazepam Glucuronide

Recovery under a simple protocol



EASY TO CLEAVE

Temazepam Glucuronide

Amitriptyline-N-β-D-Glucuronide

(+)-11-nor-9-Carboxy- D9-THC glucuronide

Oxazepam Glucuronide

HARD TO CLEAVE


Morphine-3-β-D-Glucuronide

Oxymorphone-3-β-D-Glucuronide


Codeine-6-β-D-Glucuronide

Codeine-6-β-D-Glucuronide was the hardest analyte to cleave.

2nd Under Challenging Conditions Protocol



Recovery under challenging conditions: Let's add a competitor



(+)-11-nor-9-Carboxy- D9-THC glucuronide (2,000 ng/mL)

Codeine-6-β-D-Glucuronide

500 ng/mL (open symbols)


5,000 ng/mL (solid symbols)

All samples contain 100,000 ng/mL


Acetaminophen Glucuronide as competitor

High recoveries of Codeine-6-β-D-Glucuronide are achieved even at higher analyte concentration under challenging conditions with longer incubation

3rd Realistic Conditions Protocol



Recovery under challenging conditions: higher analyte concentrations



Dihydrocodeine-6-β-D-Glucuronide

Tapentadol Glucuronide


Hydromorphone-3-β-D-Glucuronide

Morphine-3-β-D-Glucuronide

All samples contain 100,000 ng/mL

Acetaminophen Glucuronide as competitor

Recovery under challenging conditions: higher analyte concentrations



Dihydrocodeine-6-β-D-Glucuronide

Tapentadol Glucuronide


Hydromorphone-3-β-D-Glucuronide

Morphine-3-β-D-Glucuronide

All samples contain 100,000 ng/mL

Acetaminophen Glucuronide as competitor

Recovery under challenging conditions: higher analyte concentrations



Dihydrocodeine-6- β -D-Glucuronide

Tapentadol Glucuronide

Hydromorphone-3- β -D-Glucuronide

Morphine-3- β -D-Glucuronide

All samples contain 100,000 ng/mL

Acetaminophen Glucuronide as competitor

Summary Table

Type of Analyte	Concentration, ng/mL	Hydrolysis time, min	Recovery, %
Codeine	2,500	5	>85
	2,500	10	>95
	5,000	15	>95
Oxymorphone	2,500	5	91
Oxazepam	2,500	5	98
Amitriptyline	2,500	5	98
Temazepam	2,500	5	100
THC-COOH	2,500	5	100
Dihydrocodeine	20,000	15	97
Hydromorphone	20,000	15	98
Tapentadol	20,000	15	99
Morphine	20,000	15	100

What have we learned about **B-One**?

- Works for a **Comprehensive Panel of Drugs**
- Works **by cleaving** more than one drug at a time
- Works under **higher analyte** concentrations
- Works under **different extraction** methods
- Works under **different incubation** times



Take Home Message

B-One is able to **hydrolyze** conjugated drugs from **different groups** despite being challenged to non-ideal and more **real lab conditions**, demonstrating it is a good **tool** to be used for analyzing a **comprehensive panel**.

B-One has allowed labs to validate an **easier procedure** that makes sample prep simpler by **saving steps**, **saving time**, and **eliminating storage concerns** since it is stable at room temperature and refrigerated.

Finden provides a **reliable product** that is **scientifically proven** to be highly efficient.



Acknowledgements

- Nicholas Chestara, DPX Technologies, San Diego, CA.
- Jack Andrews, Kura Biotech, Puerto Varas, Chile.
- Janet Jones, Kura Biotech, Puerto Varas, Chile.
- María Jose Nuñez, Kura Biotech, Puerto Varas, Chile.
- Camila Berner, Kura Biotech, Puerto Varas, Chile.



Contact:

info@kurabiotech.com

www.kurabiotech.com

How to get a Free Sample:

To ask for a sample, please email
sales@kurabiotech.com

Kura Biotech's products are manufactured under ISO 9001:2015 International Certification and distributed to at least 15 countries.

©2021 The copyright to this document is owned by Kura Biotech S.p.A. – Kura Biotech®. No part of this document may be reproduced, stored in a retrieval system, or transmitted, in any form or by any means, electronic, mechanical, photocopying, recording, or otherwise, without their prior written permission.

KURA's products and services are protected under corresponding patents and trademark portfolio, including patent granted and pending applications, trademark registers and pending applications.

KURA
biotech